

CSF1R / MCSF Receptor / CD115 Antibody, Rabbit MAb



Catalog Number: 50059-R220

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

| | |
|---------------------|--|
| Immunogen: | Recombinant Mouse CSF1R / MCSF Receptor / CD115 protein (Catalog#50059-M08H) |
| Preparation | This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CSF1R / MCSF Receptor / CD115 (rM CSF1R / MCSF Receptor / CD115; Catalog#50059-M08H; NP_001032948.2; Met 1-Ser 511). |
| Ig Type: | Rabbit IgG |
| Clone ID: | 220 |
| Specificity: | Mouse CSF1R / MCSF Receptor / CD115 |
| | No cross-reactivity in ELISA with Human CSF1R |
| Formulation: | 0.2 µm filtered solution in PBS |
| Storage: | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |

APPLICATIONS

| | |
|----------------------|-------|
| Applications: | ELISA |
|----------------------|-------|

RECOMMENDED CONCENTRATION

| | |
|--------------|---|
| ELISA | ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CSF1R. |
|--------------|---|

Please Note: Optimal concentrations/dilutions should be determined by the end user.