

**VE-Cadherin / CD144 / CDH5
Antibody, Rabbit PAb, Antigen
Affinity Purified**

Catalog Number: 50192-RP02



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Mouse VE-Cadherin / CD144 / CDH5 protein (Catalog#50192-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse VE-Cadherin / CD144 / CDH5 (rM VE-Cadherin / CD144 / CDH5; Catalog#50192-M08H; NP_033998.2; Met 1-Ala 592). VE-Cadherin / CD144 / CDH5 specific IgG was purified by Mouse VE-Cadherin / CD144 / CDH5 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse VE-Cadherin / CD144 / CDH5
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CDH5.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.