

FCGR3 Antibody (PerCP), Rabbit MAb



Catalog Number: 50326-R002-C

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD16 / FCGR3 protein (Catalog#50326-M08H)
Reagents:	PerCP-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD16 / FCGR3 (rM CD16 / FCGR3; Catalog#50326-M08H; P08508-1; Met1-Thr215) and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.
Ig Type:	Rabbit IgG
Clone ID:	002
Specificity:	Mouse CD16 / FCGR3
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	Aqueous solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.


FCGR3 Antibody (PerCP), Rabbit MAb

Catalog Number: 50326-R002-C



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Mouse FCGR3(CD16) expression on mouse splenocytes. BABL/c splenocytes were stained with PerCP-conjugated anti-Mouse FCGR3(CD16) and FITC-conjugated anti-Mouse CD3 (BD Pharmingen™ Cat. No. 553061). The dot plots were derived from events with the forward and side light-scatter characteristics of intact cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.