



Catalog Number: 50368-R017

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Mouse MCP-1/CCL2 Protein (Catalog#50368-M08B)
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse MCP-1/CCL2 (rM MCP-1/CCL2; Catalog#50368-M08B; NP_035463.1; Met1-Asn148).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	017
<b>Specificity:</b>	Mouse MCP-1/CCL2
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	ELISA(Cap)
RECOMMENDED CONCENTRATION	
<b>Sandwich ELISA (Capture Ab)</b>	ELISA(Cap): 1:250-1:2000 This antibody will detect Mouse MCP-1/CCL2 in ELISA pair set (Catalog: # SEK50368). In a sandwich ELISA, it can be used as capture antibody.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**