

CCL2 / MCP-1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50368-T16

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CCL2 / MCP-1 Protein (Catalog#50368-M08B)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse CCL2 / MCP-1 (rM CCL2 / MCP-1; Catalog#50368-M08B; NP_035463.1; Met1-Asn148). CCL2 / MCP-1 specific IgG was purified by Mouse CCL2 / MCP-1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse CCL2 / MCP-1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:5000-1:10000
This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CCL2 / MCP-1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.