

# Clusterin / Apolipoprotein J / Apo-J / CLU Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 50485-RP02



Sino Biological  
Biological Solution Specialist

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Mouse Clusterin / Apolipoprotein J / Apo-J / CLU protein (Catalog#50485-M08H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Mouse Clusterin / Apolipoprotein J / Apo-J / CLU (rM Clusterin / Apolipoprotein J / Apo-J / CLU; Catalog#50485-M08H; NP_038520.2; Met 1-Glu 448). Clusterin / Apolipoprotein J / Apo-J / CLU specific IgG was purified by Mouse Clusterin / Apolipoprotein J / Apo-J / CLU affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Mouse Clusterin / Apolipoprotein J / Apo-J / CLU
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	ELISA
RECOMMENDED CONCENTRATION	
<b>ELISA</b>	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CLU.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**