

MMP-9 / MMP9 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50560-T16

GENERAL INFORMATION

Immunogen:	Recombinant Mouse MMP-9 / MMP9 protein (Catalog#50560-MNAH)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse MMP-9 / MMP9 (rM MMP-9 / MMP9; Catalog#50560-MNAH; P41245-1; Met1-Pro730). MMP-9 / MMP9 specific IgG was purified by Mouse MMP-9 / MMP9 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse MMP-9 / MMP9
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse MMP-9 / MMP9.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.