

CD83 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50694-RP02

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD83 protein (Catalog#50694-M08H)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse CD83 (rM CD83; Catalog#50694-M08H; O88324; Met 1-Arg 133). CD83 specific IgG was purified by Mouse CD83 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse CD83
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse CD83.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.