

Interferon Gamma / IFN gamma / IFNG Antibody, Rabbit MAb



Catalog Number: 50709-R002

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse Interferon Gamma / IFN gamma / IFNG Protein (Catalog#50709-MNAH)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse Interferon Gamma / IFN gamma / IFNG (rM Interferon Gamma / IFN gamma / IFNG; Catalog#50709-MNAH; NP_032363.1; Met1-Cys155).
Ig Type:	Rabbit IgG
Clone ID:	002
Specificity:	Mouse Interferon Gamma / IFN gamma / IFNG
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse Interferon Gamma / IFN gamma / IFNG.

Please Note: Optimal concentrations/dilutions should be determined by the end user.