

IL2/IL-2/Interleukin-2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 51061-T16

GENERAL INFORMATION

Immunogen:	Recombinant Mouse IL2/IL-2/Interleukin-2 Protein (Catalog#51061-MNAE)
Preparation	Produced in rabbits immunized with purified, recombinant Mouse IL2/IL-2/Interleukin-2 (rM IL2/IL-2/Interleukin-2; Catalog#51061-MNAE; P04351; Ala21-Gln169). IL2/IL-2/Interleukin-2 specific IgG was purified by Mouse IL2/IL-2/Interleukin-2 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse IL2/IL-2/Interleukin-2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:5000-1:10000
This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse IL2/IL-2/Interleukin-2.

Please Note: Optimal concentrations/dilutions should be determined by the end user.