

# SARS-CoV-2 Nucleocapsid Antibody, Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 68097-MM183

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant SARS-CoV-2 Nucleocapsid CTD Protein (Catalog#40588-V07E5)
<b>Preparation</b>	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant SARS-CoV-2 Nucleocapsid CTD Protein (Catalog#40588-V07E5; YP_009724397.2; Lys248-Pro365). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
<b>Clone ID:</b>	183
<b>Ig Type:</b>	Mouse IgG1
<b>Specificity</b>	SARS-CoV-2 Nucleocapsid Protein  <b>Has cross-reactivity in WB with</b> SARS-CoV-2 Nucleocapsid Protein (Cat# 40588-V08B)  <b>Has cross-reactivity in ELISA with</b> SARS-CoV-2 Nucleocapsid Protein (Cat# 40588-V07E) SARS-CoV-2 Nucleocapsid CTD Protein (Cat# 40588-V07E5) SARS-CoV-2 Nucleocapsid (aa 248-407) Protein (Cat# 40588-V07E6) SARS-CoV-2 Nucleocapsid (R203K, G204R) Protein (Cat# 40588-V07E1) SARS-CoV-2 Nucleocapsid (P13L) Protein (Cat# 40588-V07E2) SARS-CoV-2 Nucleocapsid (D3L, R203K, G204R, S235F) Protein (Cat# 40588-V07E7) SARS-CoV-2 Nucleocapsid (D3L, S235F) Protein (Cat# 40588-V07E8) SARS-CoV-2 Nucleocapsid (T205I) Protein (Cat# 40588-V07E9) SARS-CoV-2 Nucleocapsid (P80R) Protein (Cat# 40588-V07E11) SARS-CoV-2 Nucleocapsid (A220V) Protein (Cat# 40588-V07E13) SARS-CoV-2 Nucleocapsid (P199L, M234I) Protein (Cat# 40588-V07E18) SARS-CoV-2 Nucleocapsid (R203M, D377Y) Protein (Cat# 40588-V07E20) SARS-CoV-2 Nucleocapsid (A119S, R203K, G204R, M234I) Protein (Cat# 40588-V07E21) SARS-CoV-2 Nucleocapsid (P67S, P199L) Protein (Cat# 40588-V07E22) SARS-CoV-2 Nucleocapsid (P13L, T205I, A220V) Protein (Cat# 40588-V07E24) SARS-CoV-2 Nucleocapsid (P80R, A119S, R203K, G204R, M234I) Protein (Cat# 40588-V07E25) SARS-CoV-2 Nucleocapsid (P67S, P199L, M234I) Protein (Cat# 40588-V07E27) SARS-CoV-2 Nucleocapsid (D63G, R203M, D377Y) Protein (Cat# 40588-V07E29)  <b>No cross-reactivity in ELISA with</b> SARS-CoV-2 Nucleocapsid NTD Protein (Cat# 40588-V07E10)
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.
APPLICATIONS	
<b>Applications:</b>	WB,ELISA  IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
<b>Western Blot</b>	WB: 1:500-1:2000
<b>ELISA</b>	ELISA: 1:1000-1:2000

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*

Manufactured By Sino Biological Inc.


Tel: 215-583-7898 (US); +86-400-890-9989(Global) • <http://www.sinobiological.com>

# SARS-CoV-2 Nucleocapsid Antibody, Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 68097-MM183



## Anti-SARS-CoV-2 Nucleocapsid mouse monoclonal antibody at 1:1000 dilution.

Sample: SARS-CoV-2 (2019-nCoV)  
Nucleocapsid-His recombinant  
Protein(Cat#40588-V08B)  
Lane A: 10ng

## Secondary

Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.  
Performed under reducing conditions.

Manufactured By Sino Biological Inc.

Tel: 215-583-7898 (US); +86-400-890-9989(Global) ● <http://www.sinobiological.com>