

TNF-alpha / TNFA Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 70003-RP02

GENERAL INFORMATION

Immunogen:	Recombinant Canine TNF-alpha protein (Catalog#70003-DNAE)
Preparation	Produced in rabbits immunized with purified, recombinant Canine TNF-alpha (rCanine TNF-alpha; Catalog#70003-DNAE; NP_001003244.4; Val 77-Leu 233). TNF-alpha specific IgG was purified by Canine TNF-alpha affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Canine TNF-alpha
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:5000-1:10000
This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Canine TNFα.

Please Note: Optimal concentrations/dilutions should be determined by the end user.