



Catalog Number: 70014-R001

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Canine IL-2 Protein (Catalog#70014-DNAE)
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Canine IL-2 (Catalog#70014-DNAE; Q29416; Ala21-Thr155,147Cys/Ser).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	001
<b>Specificity:</b>	Canine IL-2
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	ELISA
RECOMMENDED CONCENTRATION	
<b>ELISA</b>	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Canine IL-2.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**