

IL-2 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 70014-RP01

GENERAL INFORMATION

Immunogen:	Recombinant Canine IL-2 protein (Catalog#70014-DNAE)
Preparation	Produced in rabbits immunized with purified, recombinant Canine IL-2 (Catalog#70014-DNAE; Q29416; Ala21-Thr155,147Cys/Ser). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Canine IL-2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:1000-1:2000
This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Canine IL-2.

Please Note: Optimal concentrations/dilutions should be determined by the end user.