

ACE2 / Angiotensin-Converting Enzyme 2 Antibody, Rabbit MAb

Catalog Number: 80031-R001

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Rat ACE2 protein (Catalog#80031-R08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat ACE2 (rR ACE2; Catalog#80031-R08H; Q5EGZ1; Met 1-Thr 740).
Ig Type:	Rabbit IgG
Clone ID:	001
Specificity:	Rat ACE2 Has cross-reactivity in ELISA with Human ACE2, Mouse ACE2, Rhesus ACE2.
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat ACE2.

Please Note: Optimal concentrations/dilutions should be determined by the end user.