

FAS / CD95 / APO-1 / TNFRSF6 Antibody, Rabbit MAb



Catalog Number: 80043-R105

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Rat FAS / TNFRSF6 / CD95 protein (Catalog#80043-R08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat FAS / TNFRSF6 / CD95 (rR FAS / TNFRSF6 / CD95; Catalog#80043-R08H; NP_631933.2; Met1-Lys170).
Ig Type:	Rabbit IgG
Clone ID:	105
Specificity:	Rat FAS / TNFRSF6 / CD95
	Has cross-reactivity in ELISA with Human FAS No cross-reactivity in ELISA with Mouse FAS
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat FAS.

Please Note: Optimal concentrations/dilutions should be determined by the end user.