

**MMP9/MMP-9/CLG4B Antibody,
Rabbit MAb**

Catalog Number: 80049-R213

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Rat MMP9/MMP-9/CLG4B Protein (Catalog#80049-R08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat MMP9/MMP-9/CLG4B (rR MMP9/MMP-9/CLG4B; Catalog#80049-R08H; EDL96479.1; Met1-Pro708).
Ig Type:	Rabbit IgG
Clone ID:	213
Specificity:	Rat MMP9/MMP-9/CLG4B
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat MMP9/MMP-9/CLG4B.

Please Note: Optimal concentrations/dilutions should be determined by the end user.