

MMP-9 / MMP9 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 80049-T24

GENERAL INFORMATION

Immunogen:	Recombinant Rat MMP-9 / MMP9 Protein (Catalog#80049-R08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rat MMP-9 / MMP9 (rR MMP-9 / MMP9; Catalog#80049-R08H; EDL96479.1; Met1-Pro708). MMP-9 / MMP9 specific IgG was purified by Rat MMP-9 / MMP9 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat MMP-9 / MMP9
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA, IHC-P

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat MMP-9 / MMP9.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

MMP-9 / MMP9 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 80049-T24



Immunochemical staining of rat MMP9 in rat spleen with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).