

VE-Cadherin / CD144 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 80276-T16

GENERAL INFORMATION

Immunogen:	Recombinant Rat VE-Cadherin / CD144 protein (Catalog#80276-R08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rat VE-Cadherin / CD144 (Catalog#80276-R08H; NP_001100877.1; Met1-Gln585). VE-Cadherin / CD144 specific IgG was purified by Rat VE-Cadherin / CD144 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat VE-Cadherin / CD144
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:5000-1:10000
This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat VE-Cadherin / CD144.

Please Note: Optimal concentrations/dilutions should be determined by the end user.