

**CSF1R / MCSF Receptor / CD115
Antibody, Rabbit PAb, Antigen
Affinity Purified**

Catalog Number: 80447-T16



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Rat CSF1R / MCSF Receptor / CD115 Protein (Catalog#80447-R08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rat CSF1R / MCSF Receptor / CD115 (rR CSF1R / MCSF Receptor / CD115; Catalog#80447-R08H; D4ACA7; Met1-Glu510). CSF1R / MCSF Receptor / CD115 specific IgG was purified by Rat CSF1R / MCSF Receptor / CD115 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat CSF1R / MCSF Receptor / CD115
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat CSF1R / MCSF Receptor / CD115.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.