

IFN gamma Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 90008-RP01

GENERAL INFORMATION

Immunogen:	Recombinant Rhesus IFN gamma protein (Catalog#90008-CNAE)
Preparation	Produced in rabbits immunized with purified, recombinant Rhesus IFN gamma (Catalog#90008-CNAE; NP_001028077.1; Gln24-Gln165). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rhesus IFN gamma
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Rhesus IFN gamma.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.