

# CD32b / FCGR2B / Fc gamma RIIB Antibody, Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 90014-MM06

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Cynomolgus CD32b / FCGR2B / Fc gamma RIIB Protein (Catalog#90014-C08H)
<b>Preparation</b>	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Cynomolgus CD32b / FCGR2B / Fc gamma RIIB ( Catalog#90014-C08H; AAL92097.1; Met1-Pro217). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Mouse IgG1
<b>Clone ID:</b>	06
<b>Specificity:</b>	Cynomolgus CD32b / FCGR2B / Fc gamma RIIB
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

<b>Applications:</b>	IHC-P
----------------------	-------

## RECOMMENDED CONCENTRATION

<b>IHC-P</b>	IHC-P: 1:50-1:200
--------------	-------------------

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**


# CD32b / FCGR2B / Fc gamma RIIB Antibody, Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 90014-MM06

---



Immunochemical staining of cynomolgus CD32b in cynomolgus liver with mouse monoclonal antibody (1:60, formalin-fixed paraffin embedded sections).