

**TRAIL R4 / CD264 / TNFRSF10D
Antibody, Rabbit PAb, Antigen
Affinity Purified**

Catalog Number: 90100-T24



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Rhesus TRAIL R4 / CD264 / TNFRSF10D protein (Catalog#90100-C08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rhesus TRAIL R4 / CD264 / TNFRSF10D (Catalog#90100-C08H; XP_001107922.1; Met1-Leu206). TRAIL R4 / CD264 / TNFRSF10D specific IgG was purified by Rhesus TRAIL R4 / CD264 / TNFRSF10D affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rhesus TRAIL R4 / CD264 / TNFRSF10D
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA, IHC-P
----------------------	--------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rhesus TRAIL R4 / CD264 / TNFRSF10D.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

TRAIL R4 / CD264 / TNFRSF10D
Antibody, Rabbit PAb, Antigen
Affinity Purified


Catalog Number: 90100-T24




Sino Biological
Biological Solution Specialist



Immunochemical staining of Rhesus TNFRSF10D in Rhesus liver with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



Immunochemical staining of Rhesus TNFRSF10D in Rhesus testis with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



Immunochemical staining of Rhesus TNFRSF10D in human kidney with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).