




Supreme Antibody Supplier


Leading 2nd Generation EliteRmab[®] Technology

- Fully validated in various applications
- Ideal for mouse samples
- Higher Affinity & Sensitivity than Mouse/Rat Antibodies
- High specificity & low background: perfect for IHC



Sino Biological CD71 rabbit MAb, placenta tissue IHC

VS



CD71 rabbit Polyclonal antibody (brand A), placenta tissue IHC




CD71 mouse MAb (brand B), placenta tissue IHC

Quality Assurance


- 100% produced in-house, no OEM products
- Every batch is fully validated in various applications
- Validated in Knock-out cells/tissues

IHC




Human LAMP1 in human breast carcinoma (11215-MM07)

IF




Human LAMP1 in MCF7 cells (11215-MM07)

FCM



Human LAMP1 on Jurkat cells (11215-MM07)

Human TFRC (11020-R001)



Anti-TFRC rabbit monoclonal antibody at 1:500 dilution
A: TFRC knockout HeLa Whole Cell Lysate
B: HeLa Whole Cell Lysate

Custom Antibody Development & Production Services

Rich experience: 10,000+ custom antibody projects (99% success rate)
Fast speed: Fastest delivery time in 2-3 weeks from sequence to purified antibodies
Large-scale production ability: 30-100 grams custom-tailored mAb production in weeks

- Antibody Development services
 - Rabbit Mab Development Service
 - Mouse Mab Development Service
 - Rabbit Pab Development Service
 - ELISA/WB/IHC/FACS Ab Development Service
- Antibody Production CRO Services
 - Large Scale Antibody Production Service
 - Transient Transfection Antibody Production Service
- Antibody Humanization Engineering CRO Service

Bulk Order, savings more !

- In stock & ready to ship
- Bulk order with amazing discount

How to order ?

Featured Antibody Products


1 Primary antibody

- 10,000+ antibodies, all fully validated

Research area	Antibodies
Immune Checkpoint	400+
Immunotherapy	1,100+
Targeted Therapy	1,000+
Growth Factor	200+
CD Antigens	2,500+
Enzymes	2,000+
Cytokines	800+
Receptors	2,000+
Complement System	200+

2 Reliable loading antibody


- Detect Whole-cell lysate, cytoplasmic protein, mitochondrial protein, and nucleoprotein.



3 High quality tag antibody


- Tag: His, GFP, GST, HA, Myc, DYKDDDDK (FLAG®)

His Tag Antibody (105327-MM02T):



Broad spectrum-Identify more his-tag proteins

A-P: proteins with His tag at different positions



High sensitivity-detect lower expression levels

Different loading quantity of IL2Ra-His protein sample

4 Other elite antibodies

- M13 Antibody (11973-MM05T)
- M13 Antibody (HRP) (11973-MM05T-H)

Secondary antibody

- Species: mouse, rabbit, goat
- Anti-human, mouse, rat, rabbit
- Immunogen Category: IgG, IgM, IgG (Fc)
- Conjugate: HRP, Biotin

Hot Secondary Antibodies List

Product name	Cat. No.
Goat Anti-Rabbit IgG Secondary Antibody (HRP)	SSA004
Goat Anti-Mouse IgG Secondary Antibody (HRP)	SSA007
Goat Anti-Human IgG-Fc Secondary Antibody (HRP)	SSA001
Goat Anti-Mouse IgG-Fc Secondary Antibody	SSA020
Goat Anti-Rat IgG Secondary Antibody (HRP)	SSA005

5 Leading-edge ELISA Kits

- 8 stringent QC standards: Recovery, Natural sample test, Cross-reactivity, Interference, Stability, Precision, The limit of detection (LOD), Range

Hot ELISA Kits List

Target	Sample Type	Cat. No.
Human TNF alpha	Cell culture supernates	KIT10602
Human IL-10	Cell culture supernates	KIT10947
Mouse IL-18	Mouse serum	KIT50073
Human TGF beta 1	Human serum	KIT10804
Human MCP-1/CCL2	Human serum	KIT10134
Human GM-CSF	Cell culture supernates	KIT10015
Human IL-6	Cell culture supernates	KIT10395A
Human CXCL9	Human serum	KIT10888

FLAG® is a registered trademark of Sigma-Aldrich Biotechnology.

Sino Biological US Inc.

Address: 1400 Liberty Ridge Drive, Suite 101, Wayne, PA 19087

Tel: 215-583-7898 Fax: 267-657-0217

Email: order_us@sinobiological.com

Website: www.sinobiological.com



Sino Biological
Biological Solution Specialist